



Product Datasheet

HSP70 antibody (orb767029)


Description	Rabbit polyclonal antibody to HSP70
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	Synthesized peptide derived from human HSP70 around the non-acetylation site of K246.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	70
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of human-kid...



Immunofluorescence analysis of human-kid...



Western Blot analysis of MCF7 cells usin...