


Product Datasheet

Tubulin alpha antibody (orb767035)


Description	Rabbit polyclonal antibody to Tubulin alpha
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from human Tubulin α around the non-acetylation site of K163.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	50
Expiration Date	12 months from date of receipt.



Western Blot analysis of PC12 cells using Tubulin alpha antibody (orb767035).



Immunohistochemical analysis of paraffin-embedded tissue using Tubulin alpha antibody (orb767035).



Immunohistochemical analysis of paraffin-embedded tissue using Tubulin alpha antibody (orb767035).