



Product Datasheet

IL12A p35 antibody (orb767047)


Description	Rabbit polyclonal antibody to IL12A p35
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human IL12A. AA range:31-80
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35
Expiration Date	12 months from date of receipt.



Western Blot analysis of MCF7 cells using IL12A p35 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using IL12A p35 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using IL12A p35 antibody.