


Product Datasheet

FOXP3 antibody (orb767051)


Description	Rabbit polyclonal antibody to FOXP3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the C-terminal region of human FOXP3. AA range:381-430
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	47
Expiration Date	12 months from date of receipt.



Western Blot analysis of PC12 cells using FOXP3 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using FOXP3 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using FOXP3 antibody.