



Product Datasheet

Rab 5A antibody (orb767061)


Description	Rabbit polyclonal antibody to Rab 5A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human RAB5A. AA range:71-120
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	24
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using Rab 5A antibody.



Immunohistochemical analysis of paraffin-embedded tissue using Rab 5A antibody.



Immunohistochemical analysis of paraffin-embedded tissue using Rab 5A antibody.