



Product Datasheet

Lipoprotein Lipase antibody (orb76707)

Description	Mouse monoclonal antibody to Lipoprotein Lipase
Species/Host	Mouse
Reactivity	Bovine, Gallus, Guinea pig, Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FC, IF, WB
Immunogen	Purified bovine milk lipoprotein lipase protein
Target	Lipoprotein Lipase
Preservatives	Purified antibody in PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P06858
Entrez	4023
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; IF: Use at 2 µg/ml WB: Use at an assay dependent dilution. Predicted molecular weight: 53 kDa. Flow Cyt: Use 1µg for 10 ⁶ cells.
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of using Lip...



Flow cytometric analysis of using Lipopr...