



Product Datasheet

IL1RII antibody (orb767075)


Description	Rabbit polyclonal antibody to IL1RII
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human IL1R2. AA range:141-190
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	43
Expiration Date	12 months from date of receipt.



Western Blot analysis of PC12 cells using IL1RII antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using IL1RII antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using IL1RII antibody.