


Product Datasheet

Mad1L1 antibody (orb76709)

Description	Mouse monoclonal antibody to Mad1L1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IP
Immunogen	Recombinant full length Mad1 protein (human)
Target	Mad1
Preservatives	Purified antibody containing PBS + 0.1% sodium azide
Concentration	0.9mg/ml
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q9Y6D9
Entrez	8379
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; IF: 1/500 - 1/1000. IP: Use at a concentration of 5 µg/0.5 - 1ml of extract
Expiration Date	12 months from date of receipt.



Western blot analysis of A-431 cells lys...