



Product Datasheet

TMEM173 antibody (orb767093)

Description	Rabbit polyclonal antibody to TMEM173
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from Transmembrane protein 173 at AA range: 301-350
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42,23
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa, K562, 4T1...



Immunohistochemical analysis of paraffin...