



Product Datasheet

POLR3G antibody (orb767094)


Description	Rabbit polyclonal antibody to POLR3G
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human POLR3G. AA range:151-200
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	26
Expiration Date	12 months from date of receipt.



Western Blot analysis of PC12 cells using POLR3G antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using POLR3G antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using POLR3G antibody.