


Product Datasheet

Mad1p antibody (orb76710)

Description	Rabbit monoclonal antibody to Mad1p
Species/Host	Rabbit
Reactivity	Yeast
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant fusion protein (<i>S. cerevisiae</i>)
Target	Mad1
Preservatives	Unpurified antiserum containing 0.1% sodium azide
Note	For research use only
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; WB: 1/100. Detects a band of approximately 83 kDa
Expiration Date	12 months from date of receipt.



Western blot analysis of whole cell lysa...