



Product Datasheet

TRPC3 antibody (orb767125)


Description	Rabbit polyclonal antibody to TRPC3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human TRPC3. AA range:411-460
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	97
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using TRPC3 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using TRPC3 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using TRPC3 antibody.