



Product Datasheet

Ub antibody (orb767146)


Description	Rabbit polyclonal antibody to Ub
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	Synthesized peptide derived from the N-terminal region of human Ub.
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80 50
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of human-lun...



Immunofluorescence analysis of human-lun...



Immunofluorescence analysis of mouse-spl...