



Product Datasheet

IL6Ralpha antibody (orb767151)


Description	Rabbit polyclonal antibody to IL6Ralpha
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human IL6R. AA range:221-270
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	53
Expiration Date	12 months from date of receipt.



Western Blot analysis of AD293, 4T1 cell...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...