



Product Datasheet

11 beta-HSD1 antibody (orb767165)


Description	Rabbit polyclonal antibody to 11 beta-HSD1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	Synthesized peptide derived from 11 β -HSD1 . at AA range: internal
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	35
Expiration Date	12 months from date of receipt.



Western Blot analysis of K562 cells using 11 beta-HSD1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using 11 beta-HSD1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using 11 beta-HSD1 antibody.