



Product Datasheet

IL3R alpha antibody (orb767175)


Description	Rabbit polyclonal antibody to IL3R alpha
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human IL3RA. AA range:241-290
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	42
Expiration Date	12 months from date of receipt.



Western Blot analysis of HepG2 cells using IL3R alpha antibody (orb767175).



Immunohistochemical analysis of paraffin-embedded tissue sections using IL3R alpha antibody (orb767175).



Immunohistochemical analysis of paraffin-embedded tissue sections using IL3R alpha antibody (orb767175).