

Product Datasheet

MRSA antibody (orb76719)

Description	Mouse monoclonal antibody to MRSA
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, WB
Immunogen	Native Staphylococcus aureus (methicillin resistant).
Target	MRSA (PP2A)
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	As detection antibody in ELISA we suggest a range of 1/5,000 – 1/25,000. Optimal antibody dilution should be determined by user laboratory since it may be dependent on the detection system used in both ELISA and IF
Expiration Date	12 months from date of receipt.