



Product Datasheet

CD62P antibody (orb767198)


Description	Rabbit polyclonal antibody to CD62P
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human SELP. AA range:81-130
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	91
Expiration Date	12 months from date of receipt.



Western Blot analysis of NIH-3T3, KB cel...



Western Blot analysis of 3T3 cells using...



Immunohistochemical analysis of paraffin...