


Product Datasheet

Mus81 antibody [Out of stock] (orb76720)

Description	Mouse monoclonal to Mus81.
Species/Host	Mouse
Tested Applications	IP, WB
Immunogen	Mus81 protein
Target	Mus81
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Entrez	117282
Dilution Range	WB:1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of using Mus81 ant...