



Product Datasheet

LIR-1 antibody (orb767201)

Description	Rabbit polyclonal antibody to LIR-1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human LILRB1. AA range:21-70
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	72
Expiration Date	12 months from date of receipt.



Western Blot analysis of SKOV3 cells using LIR-1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using LIR-1 antibody.