


Product Datasheet


Cadherin 5 antibody (orb767216)

Description	Rabbit polyclonal antibody to Cadherin 5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from the Internal region of human CDH5. AA range:391-440
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	86
Expiration Date	12 months from date of receipt.




Immunofluorescence analysis of A549. 1,

...



Western blot analysis of lysates from 1)...



Immunofluorescence analysis of human-lun...