

Product Datasheet

Ndfip1 antibody (orb76722)

Description	Hamster monoclonal antibody to Ndfip1
Species/Host	Hamster
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	N-terminal peptide of murine Ndfip 1 coupled to KLH (VEPACGSGYQQLQNEEPGE)
Target	Ndfip1
Preservatives	Purified antibody containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	Q8R0W6
Entrez	65113
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;WB: 1:500. Predicted molecular weight: 25 kDa
Expiration Date	12 months from date of receipt.