



Product Datasheet

Nesprin1 antibody (orb76725)

Description	Mouse monoclonal antibody to Nesprin1
Species/Host	Mouse
Reactivity	Human, Mouse, Rabbit, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant full length Nesprin 1(Human)
Target	Nesprin 1
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	Sp2/0
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q8NF91
Entrez	499010 , 23345 , 64009
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; ICC/IF: 1/200. WB: 1/400 - 1/800. Predicted molecular weight: 1011 kDa
Expiration Date	12 months from date of receipt.



Confocal immunofluorescence analysis of ...



Western blot analysis of extracts of nor...