



Product Datasheet

Tie-2 (phospho-Tyr992) antibody (orb767261)


Description	Rabbit polyclonal antibody to Tie-2 (phospho-Tyr992)
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized phospho-peptide around the phosphorylation site of human Tie-2 (phospho Tyr992)
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	125
Expiration Date	12 months from date of receipt.



Western Blot analysis of VEC cells using...



Immunohistochemical analysis of paraffin...



Immunohistochemical analysis of paraffin...