


Product Datasheet

Nesprin2 antibody (orb76727)

Description	Rabbit monoclonal antibody to Nesprin2
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IP, WB
Immunogen	GST-Nesprin 2 peptide
Target	Nesprin 2
Preservatives	Serum containing 0.1% sodium azide
Note	For research use only
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try between 1:500-1:2000
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...