


Product Datasheet

Nestin antibody (orb76729)

Description	Mouse monoclonal antibody to Nestin
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-P, IP, WB
Immunogen	A 150-amino-acid fragment from the cloned human nestin conjugated to GST
Target	Nestin
Preservatives	Purified antibody containing PBS 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P48681
Entrez	10763
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; WB: 1/5000. Predicted molecular weight: 177 kDa. ICC: 1/100 - 1/500. IHC(P): 1/500. Flow Cyt: 1/200.
Expiration Date	12 months from date of receipt.



Immunohistochemical analysis of formalin...