


Product Datasheet

Neurocan antibody (orb76730)

Description	Mouse monoclonal antibody to Neurocan
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Embryonic rat brain proteoglycans
Target	Neurocan
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	P3X myeloma
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O14594 , P55066 , P55067
Entrez	58982 , 1463
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; WB: Use at an assay dependent dilution. Predicted molecular weight: 143 kDa
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...