



Product Datasheet

SFRS14 antibody (orb767304)


Description	Rabbit polyclonal antibody to SFRS14
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from SFRS14 . at AA range: 260-340
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	120
Expiration Date	12 months from date of receipt.



Western Blot analysis of A549 cells using SFRS14 antibody.



Western blot analysis of SFRS14 Antibody using SFRS14 peptide.



Immunohistochemical analysis of paraffin-embedded tissue sections using SFRS14 antibody.