



Product Datasheet

RCL1 antibody (orb767307)


Description	Rabbit polyclonal antibody to RCL1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RCL1. AA range:324-373
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	41
Expiration Date	12 months from date of receipt.



Western Blot analysis of HUVEC cells using RCL1 antibody.



Western blot analysis of lysates from HUVEC cells using RCL1 antibody.



Western blot analysis of the lysates from HUVEC cells using RCL1 antibody.