


Product Datasheet

NHERF2 antibody (orb76732)

Description	Mouse monoclonal antibody to NHERF2
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Recombinant full length Human SLC9A3R2 fused to maltose binding protein
Target	SLC9A3R2
Preservatives	Purified antibody containing PBS + 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	NS1 myeloma
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q15599
Entrez	9351 , 116501 , 65962
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;WB: 1/100. Predicted molecular weight: 37 kDa
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of HaC...