



Product Datasheet

CNK1 antibody (orb767320)


Description	Rabbit polyclonal antibody to CNK1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CNKR1. AA range:491-540
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80
Expiration Date	12 months from date of receipt.



Western Blot analysis of COLO cells using CNK1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using CNK1 antibody.



Western blot analysis of lysates from COLO cells using CNK1 antibody.