

Product Datasheet

C/EBP Epsilon (phospho-Thr74) antibody (orb767375)

Description Rabbit polyclonal antibody to C/EBP Epsilon (phospho-Thr74)

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC-P, WB

Immunogen The antiserum was produced against synthesized peptide derived from human C/EBP-epsilon around the phosphorylation site of Thr74. AA range:40-89

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.


Note For research use only

Isotype IgG


Clonality Polyclonal

MW 34


Expiration Date 12 months from date of receipt.



Enzyme-Linked
Immunosorbent Assay
(Phosp...



Immunofluorescence
analysis of HeLa cell...



Immunohistochemistry
analysis of paraffi...