



Product Datasheet

TRIM3 antibody (orb767396)


Description	Rabbit polyclonal antibody to TRIM3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human TRIM3. AA range:1-50
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	81
Expiration Date	12 months from date of receipt.



Western Blot analysis of COLO cells using TRIM3 antibody.



Immunofluorescence analysis of NIH/3T3 cells using TRIM3 antibody.



Immunohistochemistry analysis of paraffin-embedded tissue using TRIM3 antibody.