



Product Datasheet

MALT1 antibody (orb767442)


Description	Rabbit polyclonal antibody to MALT1
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MALT1. AA range:301-350
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	90
Expiration Date	12 months from date of receipt.



Western Blot analysis of HUVEC cells using MALT1 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using MALT1 antibody.



Western blot analysis of HUVEC cell lysates using MALT1 antibody.