


Product Datasheet

Nup53p antibody (orb76745)

Description	Rabbit monoclonal antibody to Nup53p
Species/Host	Rabbit
Reactivity	Human, Rat, Yeast
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	Recombinant full length fusion protein (S. cerevisiae).
Target	Nup53
Preservatives	Unpurified antiserum containing 0.1% sodium azide
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q8NFH5
Entrez	129401
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try WB: 1/1000. Detects a band of approximately 53 kDa
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of G3 ...