



Product Datasheet

Rab 40B antibody (orb767450)


Description	Rabbit polyclonal antibody to Rab 40B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human RAB40B. AA range:229-278
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	34
Expiration Date	12 months from date of receipt.



Western Blot analysis of Jurkat cells us...



Immunohistochemical analysis of paraffin...



Western blot analysis of lysates from Ju...