



Product Datasheet

Rab 31 antibody (orb767459)

Description	Rabbit polyclonal antibody to Rab 31
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from Rab 31 . at AA range: 50-130
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	21
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using Rab 31 antibody.



Immunohistochemical analysis of paraffin-embedded tissue sections using Rab 31 antibody.