


Product Datasheet

ORC1 antibody (orb76746)

Description	Rat monoclonal antibody to ORC1
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Full length human ORC1 protein
Target	ORC1
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	P3-x63Ag.8.653
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q13415
Entrez	4998
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;WB: Use at 1:1000. Predicted molecular weight: 97 kDa
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of G3 ...