


Product Datasheet

ORC2 antibody (orb76747)

Description	Rat monoclonal antibody to ORC2
Species/Host	Rat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IP, WB
Immunogen	Full length protein (Human): ORC2 (hsORC)
Target	ORC2
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	P3-x63Ag.8.653
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q13416
Entrez	4999
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;WB: Use at 1:1000. Predicted molecular weight: 66 kDa
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts of G3 ...