



Product Datasheet

DNA pol β antibody (orb767489)


Description	Rabbit polyclonal antibody to DNA pol β
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human POLI. AA range:641-690
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	85
Expiration Date	12 months from date of receipt.



Western Blot analysis of HeLa cells using...



Western blot analysis of lysates from HeLa cells.



Western blot analysis of the lysates from HeLa cells.