



Product Datasheet

Frizzled-10 antibody (orb767490)


Description	Rabbit polyclonal antibody to Frizzled-10
Species/Host	Rabbit
Reactivity	Human, Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human FZD10. AA range:135-184
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	60
Expiration Date	12 months from date of receipt.



Western Blot analysis of COS-7 cells using Frizzled-10 antibody.



Western blot analysis of lysates from COS-7 cells using Frizzled-10 antibody.



Immunohistochemical analysis of paraffin-embedded tissue using Frizzled-10 antibody.