



Product Datasheet

MAN1B1 antibody (orb767499)

Description	Rabbit polyclonal antibody to MAN1B1
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Synthesized peptide derived from MAN1B1 . at AA range: 100-180
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	80
Expiration Date	12 months from date of receipt.



Western Blot analysis of 3T3 cells using...



Western blot analysis of MAN1B1 Antibody...