

Product Datasheet

Galactosidase beta 1L3 antibody (orb767508)

Description Rabbit polyclonal antibody to Galactosidase beta 1L3

Species/Host Rabbit

Reactivity Human

Conjugation Unconjugated

Tested Applications ELISA, IHC-P

Immunogen The antiserum was produced against synthesized peptide derived from human BGLL3. AA range:291-340

Preservatives PBS with 0.02% sodium azide and 50% glycerol pH 7.4.

Concentration 1 mg/ml


Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Note For research use only

Isotype IgG

Clonality Polyclonal

Expiration Date 12 months from date of receipt.



Immunohistochemistry analysis of paraffi...