



Product Datasheet

A Cyclase IX antibody (orb767526)


Description	Rabbit polyclonal antibody to A Cyclase IX
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC-P, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human ADCY9. AA range:137-186
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	150
Expiration Date	12 months from date of receipt.



Western Blot analysis of HuvEc cells usi...



Western blot analysis of ADCY9 Antibody....



Immunohistochemistry analysis of paraff...