

Product Datasheet

OVCA2 antibody (orb76753)

Description	Rabbit monoclonal antibody to OVCA2
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide: GALRKALRGRAELVC conjugated to KLH, corresponding to amino acids 27-41 of Human OVCA2
Target	OVCA2
Preservatives	Antiserum containing 0.1% sodium azide
Note	For research use only
Clonality	Monoclonal
Uniprot ID	Q8WZ82
Entrez	124641
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;WB: 1/1000. Predicted molecular weight: 24 kDa
Expiration Date	12 months from date of receipt.