


Product Datasheet

p22phox antibody (orb76754)

Description	Mouse monoclonal antibody to p22phox
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, WB
Immunogen	Partially purified human neutrophil flavocytochrome b (heparin Ultrogellectin affinity purification in octyl glucoside)
Target	Cytochrome b245 Light Chain
Preservatives	Affinity purified antiserum containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	P3U1 myeloma
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	p13498
Entrez	1535
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try; 1µg/ml for WB, 1:100 for IF, 1:100 for Flow Cytometry
Expiration Date	12 months from date of receipt.



Western blot analysis of human spleen ti...