



Product Datasheet

MRGX4 antibody (orb767543)


Description	Rabbit polyclonal antibody to MRGX4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human MRGX4. AA range:271-320
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	33
Expiration Date	12 months from date of receipt.



Western Blot analysis of HUVEC cells using MRGX4 antibody.



Immunofluorescence analysis of A549 cells using MRGX4 antibody.



Western blot analysis of lysates from HUVEC, COLO 205, and COLO 205 cells using MRGX4 antibody.