


Product Datasheet

p22phox antibody (orb76755)

Description	Mouse monoclonal antibody to p22phox
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Tissue / cell preparation (Human): Partially Purified human neutrophil flavocytochrome b (heparin Ultrogel/lectin affinity purification " in Triton X-100)
Target	Cytochrome b245 Light Chain
Preservatives	Purified antibody containing PBS 0.1% sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q01062
Entrez	1535 , 21577
Dilution Range	Optimal antibody dilution should be determined by titration, however as a guideline try;WB: Use at 1:1000. Predicted molecular weight: 21 kDa
Expiration Date	12 months from date of receipt.



Western blot analysis of human spleen ti...