



Product Datasheet

CLN5 antibody (orb767557)


Description	Rabbit polyclonal antibody to CLN5
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	The antiserum was produced against synthesized peptide derived from human CLN5. AA range:171-220
Preservatives	PBS with 0.02% sodium azide and 50% glycerol pH 7.4.
Concentration	1 mg/ml
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48
Expiration Date	12 months from date of receipt.



Western Blot analysis of 293T cells using CLN5 antibody.



Western blot analysis of CLN5 Antibody. ...



Western blot analysis of the lysates from cells expressing CLN5.