


Product Datasheet

PABP antibody (orb76757)

Description	Mouse monoclonal antibody to PABP
Species/Host	Mouse
Reactivity	Gallus, Hamster, Human, Rabbit, Xenopus
Conjugation	Unconjugated
Tested Applications	ELISA, FC, ICC, IP, WB
Immunogen	Recombinant PABP (Human) expressed from its 1.85 kbp cDNA, NcoI to SspI
Target	PABP
Preservatives	Purified antibody (from supernatant) containing PBS + 0.1% sodium azide;
Concentration	1mg/ml
Note	For research use only
Application notes	SP2/0
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	P11940
Entrez	26986
Dilution Range	Optimal antibody dilution should be determined by titration, however try:WB: Use at an assay dependent dilution. Detects a band of approximately 70 kDa (predicted molecular weight: 71 kDa)ICC/IF: Use at an assay dependent dilution. PubMed: 17977970Flow Cyt: Use 1µg for 10 ⁶ cells
Expiration Date	12 months from date of receipt.



Confocal
immunofluorescence
analysis of ...